



06701-9383-0  
info@ibs-scherer.de

## Overview

# IBS cleaning booster additive for WAS 10.100 & WAS 10.500, 5 litres

Product number: 2050426





## Description

### Description

Wherever WAS 10.100 and WAS 10.500 are already in use. In a few cases, an additional cleaning booster may be necessary.

### Advantages

- Liquid concentrate on an aqueous basis - no powder
- Gentle on materials - suitable for sensitive surfaces / workpieces (painted surfaces, brass, copper, zinc, nickel, chrome, aluminium, plastics of all kinds)
- Environmentally friendly and user-friendly
- silicate-free

### Recommended dosage

Concentration: 0.3-0.5% (based on the bath volume). For example, if your automatic washing machine contains 110 litres (IBS-MAXI 78), you will need to add 0.33 - 0.55 litres of IBS cleaning booster additive.